

POCKETPRO H²

HANDHELD ULTRASOUND



500L Specification

Frequency Range: 6-12 MHz

Maximum Depth: ~ 6cm

Weight: 260 grams (9 oz)

Height: 7 Inches

Supported Modes: B-mode, CF, PW, M

Run Time: 1hr 30 min (Scan time)
12 hours (Stand-By)

Application: MSK, Vascular Breast, Thoracic/Plural motion, Fluid Detection, Needle Guidance

Battery: Includes 2 Batteries with AC Charger
(Use one and charge the other)
Run time 1.5 hours
Charge time 1.5 hours

Image Storage: Send images to PACS: DICOM / Jpeg

Portability, Efficiency, Safety: Multiple Offices, Nursing Homes, House Calls, Sports Events

MSK Ultrasound: Scan Muscle, Cortical Bone, Tendon, Ligaments, Soft Tissue Masses, Nerve Bursa, Fascia, Cysts, and more

OS SUPPORT

Android 7.0 and above
iOS 12 and above
WINDOWS-10 (64-Bit)




Android (Samsung)

SAMSUNG Galaxy S6
SAMSUNG Galaxy Note5, FE
SAMSUNG Galaxy A5, A7 2016
SAMSUNG Galaxy J3, J5, J7 2017
SAMSUNG Galaxy Tab S3 (2017)
SAMSUNG Galaxy Tab 10.1 2016

iOS (Apple)

iPhone 6s iPad Air 2
iPad 5 iPad Pro 1
iPad Mini 6

POCKETPRO H² HANDHELD ULTRASOUND



2020IMAGING.NET/POCKETPRO