

*Universal*

*AVChoice*

Veterinary Systems



# System Overview

## Analog Systems Ready for DR

### The Veterinarian's Choice

Universal AVChoice systems deliver comprehensive veterinary imaging solutions effectively designed for efficiency while maintaining affordability.

### Built Strong and Ready

The sturdy welded construction of the table can support patients up to 400 lbs. and feature smooth surfaces with urine trap perimeters.


### Compact and Clean

The AVChoice table cabinet can be used to house system generator components greatly reduces floor space needed for installation. Complete system footprint is just 39" x 56" (with stationary or 2-way float top table option).

### Hands-Free Operation

Conveniently placed exposure switch and table top release strips are foot activated to ensure safety for the veterinarian and animals being examined. Release strips are located at both ends and table front to maximize operability from all sides of the table.

Table lock release functions are electrically powered and quiet to prevent animals from being startled during examinations.



#### STANDARD AVCHOICE SYSTEM

##### TABLE

- Stationary Table Top

##### TUBE STAND

- Fixed SID, 40" to image receptor
- Articulating arm on side of column for mounting generator control console

##### GENERATOR (single or three phase)

- 30 kW, high frequency, single phase, 400 mA, 125 kVp, with anatomical programming console
- 32 kW, high frequency, three phase, 500 mA, 125 kVp, with anatomical programming console

##### X-RAY TUBE

- 140 kHu, 125 kVp, 1.0-2.0 mm focal spots

#### SYSTEM OPTIONS

AVChoice systems offer a variety of different component options to meet the needs of any veterinary practice

##### TABLE

- 2-way float top
- 4-way float top

##### TUBE STAND


- Variable SID, 40" vertical range with electric locks

##### GENERATOR

- 40 kW and 50 kW generator options available upon request
- Pedestal stand or wall mount for control console


##### X-RAY TUBE

- 200 kHu, 125 kVp, 0.6-1.5 mm focal spots




# Technical Specs


Front View




Side View (Fixed SID)




Side View (Variable SID)



Stationary and Two-way Float Top



Four-way Float Top



[www.delmedical.com](http://www.delmedical.com)

**New York**  
28 Calvert Street  
Harrison, NY 10528  
(Tel) 800.261.9808  
(Fax) 914.835.6111

**Illinois**  
241 Covington Drive  
Bloomington, IL 60108  
(Tel) 800.800.6006  
(Fax) 847.288.7011

**DEL MEDICAL**

#9000-AVCHOICEVS 02.19