



RadmediX 2510 Lance Drive, Kettering, OH 45409 USA
Tel: 844.723.6334 | Fax: 513.278.0150 | info@radmedix.com
© 2020 RadmediX Inc. All rights reserved.



## FIRST CHIROPRACTIC X-RAY ADVANCEMENT IN DECADES

## Acuity CDR, Advanced X-Ray Technology

Radmedix has expanded our proven straight arm line with a new advanced x-ray system designed specifically for chiropractic care, the Acuity CDR. The Acuity CDR is a smaller more compact system which comes with a variety of features such as full software integration, magnetic lock brakes, 40kw generator, LED collimator, and a removeable grid.

This x-ray solution can fit in virtually any x-ray room as well as drastically cut down installation time reducing clinical downtime.



# **ACUITY CDR KEY FEATURES**

### **BUCKY ENCLOSURE**



Features powerful magnetic lock brakes from handle control to raise and lower the system effortlessly. Fully counterbalanced, ultra small footprint solution.





#### **TUBE HEAD**



Features a dual laser rotating LED collimator with precise alignment and a fixed 40in. SID system with optional Automatic Exposure Control (AEC).

## **HEAD TO TOE POSITONING**



The Acuity CDR can seamlessly raise and lower to accommodate all standing exams from skull to lower extremities. This maintains constant tube and bucky alignment, improving workflow.

